

**It's time to switch "Cides"!**



**Opti-Cide<sup>3</sup>**

**CALL FOR  
A FREE SAMPLE**

**NEW EPA APPROVALS !**

**Disinfects clinical equipment surfaces in 2 minutes!**

**Sanitizes non-clinical surfaces in 10-seconds!**

**Kills Bacteria, Viruses, Mycobacteria TB and Pathogenic Fungus.**

**Cleans dried blood, particulate matter, spores, cellular debris and other bodily soils from surfaces quickly and safely!**

**Approved for use with both reusable and single-use Microfiber cloths and mops**

- Use on infusion pumps, respiratory care equipment, OR and ER surfaces and other patient-care equipment as well as telephones, computer key boards, etc.
- Safe for clear plastics, vinyl, Plexiglas®, all metals and painted surfaces (One product for all surfaces!)
- Ready-to-Use for disinfection or sanitization. (No mixing - No mistakes!)
- Concentrated for floor cleaning and leaves no film. (Economical!)
- Available in “no-airborne mist” spray bottles for easy application (Safe to use!)
- Spray bottles may be refilled from several available larger sizes (Saves money!)
- Also available on large pre-wet single use wipes (100 per canister)



Exceeds EPA, CDC and OSHA requirements.  
Samples and Technical data are available on request.

**Call (toll free) 1-888-253-2536**



*Where innovation  
and knowledge  
are our most  
important products.*